



## Swedish Hospital Supply AB

### P-BAG Product Datasheet

<b>Issuer:</b>	<b>Date:</b>	<b>Edition:</b>	<b>BASIC UDI-DI:</b>	<b>Designation:</b>
Swedish Hospital Supply AB, (SHS AB)	2022-04-07.	2.	B-17350026250188	P-BAG
<b>Production:</b>	<b>Responsible Production:</b>	<b>Quality</b>	<b>Product holders:</b>	<b>Supplier:</b>
Fristad Plast AB (Fp)	Nicklas Forsström.	Nathaniel Forsström,	SHS AB	SHS AB

### Product description

Urine collection bag with funnel (UNISEX).

### Intended use

**P-BAG** is a very good and convenient alternative to other urine collecting aids in Hospital wards, patient transports, for wheelchair users and in home care. The P-BAG comes with a removable and reusable UNISEX funnel which can be used in most situations and positions (lying down, sitting, standing).

P-BAG collects a maximum of approximately 1 500 ml of urine and is recommended to be used during no more than 24h.

### Characteristics and benefits

- P-BAG is a urine collecting bag in a closed environment that can replace other collection vessels and suits both men and women in most situations (lying, sitting, standing).
- The P-BAG has a carrying handle that makes it easy to grip and easy to hang up.  
P-BAG can be used several times for a maximum of 24 hours thanks to a unique non-return valve.
- P-BAG is a convenient aid for the user and working-aid for the staff. Prevents direct contact with urine and odors are retained in the bag. P-BAG is made in Sweden, and from the material Polyethylene.
- P-BAG saves time and money as it can be used several times a day in most situations and is easy to manage for the patient.



## Swedish Hospital Supply AB

### Bag

Material: Polyethylene foil, LDPE, transparent

Length: Outer edge Weld. 320mm + 5mm

Width: Outer edge Weld. 183mm + 5mm

Thickness: 0.10mm + 0.01mm

Print: According to the picture or attached sample.

Color: Green

### Check valve

Material: Speedflex EMS No. 3-05. According to the attached sample.  
Punched in – top foil printed part – of the bag that interacts with injection molded funnel receiver, welded on the bag.

### Funnel receiver

Material: Injection moulded polyethylene see picture or attached sample.

### Funnel

Material: Injection moulded polyethylene see attached sample.

Handle direction: Recommended handle straight up, see picture.  
Adjustable for patient fit.

### Packaging

Article nr	UDI-DI-Carton	Number/carton	Weight/carton
100005-A	17350026250065	30 Bags/30 Funnels	approx. 1.6 kg
100006-A	17350026250072	100 Bags/10 Funnels	approx. 2,5 kg
100009-A	17350026250218	100 Funnels	approx. 4 kg
100010-A	17350026250096	150 Bags	approx. 1.5 kg
100012-A	17350026250089	10 Bags/1 Funnel	approx. 250 g



## Swedish Hospital Supply AB

100012-A	17350026250133	10 Packages/Carton	approx. 2,5 kg
----------	----------------	--------------------	----------------

### Package

Material: Unbleached corrugated carton with inner bag of transparent polyethylene.

#### **10 0005-A**

Outer dimensions length: 380 mm  
Outer dimensions width: 270 mm  
Outer dimensions height: 250 mm  
Number of details/cartons: 30 Bags/30 Funnels  
Tot. Number of cartons/pallet: 56 pcs.

#### **10 0006-A**

Outer dimensions length: 400 mm  
Outer dimensions width: 300 mm  
Outer dimensions height: 400 mm  
Number of details/cartons: 100 Bags/10  
Tot. Number of cartons/pallets: 40 pcs.

#### **10 0009-A**

Outer dimensions length: 380 mm  
Outer dimensions width: 270 mm  
Outer dimensions height: 250 mm  
Number of details/cartons: 100  
Tot. Number of cartons/pallet: 56 pcs.

#### **10 0010-A**

Outer dimensions length: 380 mm  
Outer dimensions width: 270 mm  
Outer dimensions height: 250 mm  
Number of details/cartons: 150 bags  
Tot. Number of cartons/pallet: 56 pcs.

#### **10 0012-A**

Outer dimensions length: 400 mm  
Outer dimensions width: 300 mm  
Outer dimensions height: 400 mm  
Number of details/carton: 100 Bags/10 Funnels, 10 packages with 10 Bags/ 1 Funnel.  
Tot. Number of cartons/pallet: 40 pcs.

Carton labelling: UDI-DI for each packaging volume, LOT, date of manufacture (year, month, day), expiration date (year, month, day).



## Swedish Hospital Supply AB

### **Production**

Procedures: The product is produced in accordance with ISO 9001.

### **Environmental**

Procedures: The product is produced in accordance with ISO 14001.

### **Country of origin**

Sweden

### **Sterile product**

No

### **Shelf life**

4 years, unbroken packaging

### **Storage**

Store dry at room temperature.  
Avoid exposure to high temperatures above 55C.

### **Handling Waste**

According to local regulations. Plastic recycling or incineration.

### **Other information**

HS Code/Commodity code: 39269097

<https://swehospitalssupply.se/en/>



Swedish Hospital Supply AB



Secure-Innovative-Sympathetic