



Swedish Hospital Supply AB

P-BAG SAP Datasheet

Issuer:	Date:	Edition:	BASIC UDI-DI:	Designation:
SHS, Swedish Hospital Supply, AB	2022-04-07.	2.	B-17350026250195	P-BAG-SAP
Production:	Responsible Production:	Quality Manager:	Product holders:	Supplier:
Fristad Plast AB (Fp)	Nicklas Forsstrom.	Nathaniel Forsstrom.	SHS AB	SHS AB

Product description

Urine collection bag with funnel (UNISEX) and super absorbent plate in sealed bag.

Intended Use

Urine collection bag "P-BAG" with Superabsorbent, is an extraordinarily good alternative to other urine collecting aids in situations like e.g., transports, emergency care for wheelchair users, hospitals and in home care. The P-BAG SAP comes with a detachable UNISEX funnel and is used in most situations and positions (lying down, sitting, standing).

Characteristics and benefits

- **"P-BAG" with Superabsorbent**, has a carrying handle that makes it easy to grip and to hang up. It can be used several times a day thanks to a unique check-valve and an "absorbent plate" which is reduced in size in relation to the amount of incoming urine = Rapid absorption up to about 1.2 – 1.5 liters.
- **"P-BAG" with Superabsorbent**, makes it easier for staff to avoid direct contact with urine and reduces odors. The bag and funnel are recyclable, and the absorbent managed in the same way that healthcare manage diaper waste.
- **"P-BAG" with Superabsorbent**, can transport urine containing – environmentally harmful substances such as pharmaceutical residues in absorbed form to destruction.
- **"P-BAG" with Superabsorbent**, is made in Sweden from Polyethylene.
- **"P-BAG" with Superabsorbent**, saves time for staff and money as it can be used several times in a day or at a maximum volume of approximately 1500 ml.



Swedish Hospital Supply AB

Bag

Material: Polyethylene foil, LDPE, green print

Length: Outer edge Weld. 320mm + 5mm

Width: Outer edge Weld. 183mm + 5mm

Thickness: 0.10mm + 0.01mm.

Print: Green

Check valve

Material: Speedflex EMS No. 3-05. According to the attached sample.
Punched in – green printed over foil – of the bag that interacts with injection molded funnel-receiver, welded on the bag.

Funnel receiver

Material: Injection moulded polyethylene see picture or attached sample.

Funnel

Material: Injection moulded polyethylene see attached sample.

Handle direction: Recommended, straight up. Adjustable for patient fit.

Superabsorbent

Material: About 28g ± 3g starch polymer fixed between 2 layers virgin cellulose.

Length; 250 mm,

Width; 100 mm,

Thickness; 2 mm.



Swedish Hospital Supply AB

Packaging

Article No.	UDI-DI-Cardboard	Number/carton	Weight/carton
200000-B	17350026250171	20 Bags/2 Funnels	about 1.5 kg
200001-B	17350026250126	20 Bags/20 Funnels	about 1.5 kg
200002-B	17350026250119	75 Bags/75 Funnels	about 4.5 kg

Material Cardboard: Unbleached corrugated board with a transparent polyethylene inner bag.

Article No 200001-B / 200000-B

Outer dimensions length: 380 mm.

Outer dimensions width: 270 mm.

Outer dimensions height: 250 mm.

Number of details/carton: 20 pcs.

Number of cartons/pallet layers: 8 Pcs.

Number of layers/pallet: 7 pcs.

Tot. Number of carton/pallets: 56 Pcs.

Article No 200002-B

Outer dimensions length: 400 mm.

Outer dimensions width: 300 mm.

Outer dimensions height: 400 mm.

Number of details/carton: 75 Pcs.

Number of cartons/pallet layer: 8 Pcs.

Number of layers/pallets: 5 pcs.

Tot. Number of cartons/pallets:40 Pcs.



Swedish Hospital Supply AB

Carton labelling: UDI-DI for each packing volume, LOT, date of manufacture (year, month, day), expiration date (year, month, day).

Production

The product is produced in accordance with ISO 9001.

Environment

The product is produced in accordance with ISO 14001.

Country of origin

Sweden

Sterile Product

No

Shelf life

4 years, unbroken packaging

Storage

Store dry at room temperature.
Avoid exposure to high temperatures above 55C

Waste Management

The absorbent material in the bag must not be flushed down the toilet or sewer, should be sorted as waste such as diapers.
Plastic parts can be sorted according to local regulations, plastic recycling, or incineration.

Other information

HS Code/Commodity code: 39269097

<https://swehospitalsupply.se/en/>



Swedish Hospital Supply AB



Secure-Innovative-Sympathetic