



Swedish Hospital Supply AB

## Combi-Spoon – Product Data Sheet

Issuer:	Date:	Edition:	BASIC UDI-DI:	Designation:
SHS, Swedish Hospital Supply, AB	2023-10-02	1.	17350026250300EK	Combi-Spoon

### Product description

Combi-Spoon is a cutlery with a spoon, fork and knife in one. All cutlery gathered in one can help patients/users to independently eat their meals despite impaired function.

Combi-Spoon helps patients/users and healthcare professionals to meet some of the problems that may arise at mealtimes.

Combi-Spoon is also a good and flexible cutlery to have with you in outdoor and camping situations.

### Intended Use

Combi-Spoon helps patients/users who have impaired function in their hands or arms to gather all the cutlery to one hand. In this way, the patient/user can easily divide, collect and capture the food with the Combi-Spoon.

The simple design and combined function are also appreciated in outdoor and camping situations.

### Features and Benefits

The Combi-Spoon has a combined fork and spoon-shaped top that transitions into slightly sharper edges on the sides. The edges easily divide most of the commonly used ingredients that go into a meal. The food can then be collected with the fork or spooned up with the Combi-Spoons top.



Secure-Innovative-Sympathetic



## Swedish Hospital Supply AB

### **Material**

Material: ABS-plastic  
Manufacturing: Injection Molded

### **Packaging**

Material cardboard: Unbleached corrugated cardboard.  
Labeling carton: UDI-DI for each packaging volume, LOT, date of manufacture (year, month, day), expiry date (year, month, day).

### **Production**

The product is produced in accordance with the policies and quality management systems adopted by Swedish Hospital Supply AB and ISO 90001. Quality testing is routinely carried out.

### **Country of origin**

Sweden

### **Sterile product**

No

### **Shelf life**

10 years, unbroken package

### **Storage**

Store dry, room temperature.

### **Waste management**

According to local regulation. Recycle hard plastic product.

### **Other information**

HS Code/Commodity code: 39269097

<https://swehospitalsupply.se/en/>